



## Chemical Fume Hood

*Chemical Laboratory fume hoods are used to protect the laboratory environment and the operator while working with chemicals. The fume hood actively protects the operator from inhaling toxic fumes . By installing a suitable filter it can also protect the environment.*

### Technical Specification:

<b>Width:</b>	100 , 130 , 160cm
<b>Material:</b>	Body: electro statically coated metal / Tray: SS-304L
<b>Exhaust Fan flow:</b>	1620 m <sup>3</sup> /hr - 954 cfm
<b>Exhaust Filter(Optional):</b>	Nano Carbon
<b>Class/Type:</b>	Full Exhaust
<b>Additional Options:</b>	SS Sink and water valve, gas line, second exhaust duct-fan , inside power outlet