



## CO<sub>2</sub> Incubator **INCUBE-120 CO<sub>2</sub>**

*One of the most important characteristics of CO<sub>2</sub> incubators is the regulation of moisture and carbon dioxide during cell culture. by using this device, the temperature and pH of the media can be adjusted for cell culture. Carbon dioxide incubator is used for microbiological, biochemical, genetic, pharmaceutical and biology researches. This economic CO<sub>2</sub> incubator is designed for research and development in pharmaceutical companies and academic centers. Co<sub>2</sub> and temperature stability in this device makes it a good choice for IVF and artificial insemination.*

### Technical Specification:

Working volume:	120L
Material:	Body: electro statically coated metal / Chamber: SS-304L electro-polished Ra ≤0.6µm
Main controller:	TDE PLC system
Display:	3.5" blue backlight characteric LCD
Temperature Range:	RT+3 _ 50°C
Humidity Range:	0 _ 100%RH
CO <sub>2</sub> Range:	0 _ 20%

No. 4, East Second Alley, North Mottahari Bld.,  
Safadasht, TEHRAN, IRAN  
[www.kavoshteb.com](http://www.kavoshteb.com) [info@kavoshteb.com](mailto:info@kavoshteb.com)  
tel: +98 21 65536644

**KAVOSH TEB**