



Standard Incubator

INCUBE Series

Standard incubators are suitable for a variety of routine applications where precise and consistent temperature control is required.

Also the **INCUBE-750 R** has refrigerated technology which offers precise temperature control for reliable results in pharmaceutical, food, cosmetic and microbiology research.

Technical Specification:

Working volume:	120 , 750L
Material:	Body: electro statically coated metal / Chamber: SS-304L electro-polished Ra \leq 0.6 μ m
Main controller:	TDE PLC system
Display:	3.5" blue backlight characteric LCD
Temperature Range:	RT+3 _ 50°C 10 _ 50°C for refrigerate model: 750R
Humidity Display:	0 _ 100%RH



No. 4, East Second Alley, North Mottahari Blv.,
Safadasht, TEHRAN, IRAN
www.kavoshteb.com info@kavoshteb.com
tel: +98 21 65536644

KAVOSH TEB

